

VT-1980

Electronic personal scale

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ENGLISH

EXCESS WEIGHT SCALES ANALYZER

The Bioelectric Impedance Analysis "BIA" implemented in the appliance allows calculating fat, water, muscle mass and bone mass percentage of the body with high accuracy taking into account personal data (age, gender, height) stored in memory.

DESCRIPTION

1. Display
2. Upper contact pads
3. UNIT – measuring units selection button (kg/st/lb)
4. Lower contact pads
5. ▼ - button for decreasing the displayed figure
6. The "ON" button
7. ▲ - button for increasing the displayed figure
8. SET – selection and confirmation of the data entered
9. Supporting legs
10. Battery compartment lid

SAFETY MEASURES

Before using the unit read these instructions carefully.

- The use of excess weight analyzer is **ABSOLUTELY** prohibited for pregnant women.
- The use of excess weight analyzer is **ABSOLUTELY** prohibited for persons having implanted cardio stimulators, as the operation principle of the appliance is based on running of electric current through the body of a person, it may cause malfunction of the cardio stimulator.
- Place the scales on the flat, dry and slip-proof surface.
- In order to avoid inaccuracy while weighing, do not use the scales on carpets and carpeted floor.
- To get maximally precise results while weighing put the scales each time on the same place.

- Do not stand on the scales with wet feet or on the wet surface of the scales. You may slip down and get injured.
- Be careful while standing on the scales: firstly put your one foot on the scales, make sure that the scales are stable and then put the other foot on the scales.
- Stand on the scales and distribute your body weight evenly. Stay steel while weighing, do not jump on the surface of the scales.
- Use the scales away from heating devices.
- Do not use and do not keep the scales in rooms with high humidity (more than 80%), protect the body of the scales from water or other liquids, it may influence the precision of weighing and lead to outrage.
- Be especially careful if the unit is used by children or disabled persons.
- Do not allow children to use the unit as a toy.
- If you are not going to use the scales for a long time, you should take the batteries out of the battery compartment.
- While inserting batteries, strictly follow the polarity according to the marks.
- Do not use different batteries together (new and old ones, different types of batteries etc.);
- Replace the batteries in time.
- Do not dismantle and do not repair the scales yourself. If there are any damages in the unit, apply to the authorized service center.
- The unit is intended for household usage only and is not intended for medical or commercial usage.

RECOMMENDATIONS

- Perform weighing at the same time of the day.
- Percentage of water should be weighed in the evening before supper. 2-3 hours should pass after the last meal.

- 10-12 hours should pass after the intensive physical exercises.

Attention!

- **Data concerning percentage of fat, water and muscles can be obtained only if you stand on the scales bare-foot.**
- Percentage of fat depends on different factors and may vary during a day.
- When determining the percentage of fat and other parameters of children younger than 10 years old the results may be inaccurate.
- Professional sportsmen, persons, liable to edema, persons with non-standard body proportions may get inaccurate results of percentage of fat, muscles and bone mass.
- The scales can be used for ordinary weighing of children of any age whose weight is at least 2 kg.

OPERATION OF THE SCALES

If batteries are supplied separately and when replacing the batteries:

- Open the battery compartment lid (10), by pressing the clamp.
- Insert 4 “AAA” (4x1,5V) batteries into the battery compartment, strictly following the polarity.
- Close the battery compartment lid (10).
- If you need to erase all data from the memory of the appliance, open the lid (10), remove one battery for some seconds, after it put it on its place and close the battery compartment lid (10).

If batteries are installed in the battery compartment:

- Open the battery compartment lid (10), remove the insulating pad, and close the lid (10).

MEASURING UNITS SELECTION

To select the measuring units, switch on the scales, by pressing the (6) button, se-

lect measuring units, touching the UNIT (7) button:

- “kg” - weight will be indicated in kilograms;
- “st” - weight will be indicated in stones;
- “lb” - weight will be indicated in pounds;

ATTENTION!

Meaning of the symbols on the display (1):

- In case of overloading the scales the “Err” symbol will appear on the display (1), get off the scales immediately to prevent the damage (pic.1).
- If “Lo” appears on the display, it is necessary to replace the batteries (pic.2) (see the OPERATION OF THE SCALES section).
- If during calculation of percentage of fat and water “Err” appears on the display (1) it means that your feet do not touch the upper (2) or lower (4) contact pads (pic.3).
- Shall the “ErL” symbol appear, it means very small percentage of fat in the tissues of the body that is why it is impossible to perform calculations (pic.4). Recommendations: you should stick to more nutritious dietary regime and be more attentive to your health.
- “ErH” means excess of percentage of fat in the tissues of the body that is why it is impossible to perform calculations (pic.5). Recommendations: you should stick to low-calorie dietary regime, go in for sports and be more attentive to your health.

ORDINARY WEIGHING

1. Switch on the unit by pressing the (6) button.
2. Wait until “0.0 kg”, “0:0.0 st” or “0.0 lb” appears on the display (1), depending on the unit selected by the UNIT (3) button (pic.6).
3. Stand on the scales, with your feet parallel and evenly distribute the weight.

ENGLISH

Stay motionless while weighing. Your weight will be displayed (1) during 15 seconds, after that the scales will switch off automatically (pic.7)

PERSONAL DATA INPUT

(as per the example of indicating the weight in "kg", the height data input in "cm")

1. Switch on the unit by pressing the (6) button.
2. Select the measuring unit by touching the UNIT (7) button:
"kg" - weight will be indicated in kilograms and the height in centimeters
"st" - weight will be indicated in stones and the height in feet and inches
"lb" - weight will be indicated in pounds and the height in feet and inches.
3. Touch the SET (8) button, number of the user "P0" will appear on the display (pic.8).
4. Touch the ▼ (5) or ▲ (7) button to select the user number. You may enter data for 10 users - 0,1, ...9. Confirm the selection of the user by touching the SET (8) button.
5. Touch the ▼ (5) or ▲ (7) button to select the gender symbol: ▼ - for men, ▲ - for women, and ⚡, if you go in for sports (pic.9, 10). Confirm the selection of the gender by touching the SET (8) button.
6. Then the flashing "age" symbol will appear on the display (1), touch the ▼ (5) or ▲ (7) button to enter your age in the range from 10 till 80 years (pic.11). Confirm the selection of the age by touching the SET (8) button.
7. Touch the ▼ (5) or ▲ (7) button to enter your height from 100 cm till 240 cm (pic.12). Confirm the selection of the height by touching the SET (8) button.
8. After entering all the data "0.0 kg" will appear on the display and the scales will switch to the excess weight measuring mode (pic.13).

Note:

- If time range between each personal data input is more than 10 sec, the appliance will switch off automatically.
- If no data is entered in the memory the appliance will measure fat, water, muscles and bone mass percentage of a 25 years old male 1.65 m in height.

CALCULATIONS ACCORDING TO THE WEIGHING RESULTS

1. Touch the SET (8) button to switch on the scales.
2. Touch the ▼ (5) or ▲ (7) button to select the user number and confirm it by touching the SET (8) button.
3. Subsequently touching the SET (8) button you can view the personal data (gender, age, height).
4. After indication of the height data, "0.0 kg" and the user number "P0" will appear on the display (pic.13).
5. Stand on the scales barefoot in such a way that your feet and heels are on the upper and lower contact pads (2 and 4), evenly distribute your weight (pic.14). Stay motionless while weighing. Your weight will appear on the display (1), when the moving symbol "o" appears, calculation of the percentage of fat "FAT", water "TBW", muscles mass "MUS" and bone mass "BONE" in the human body is being performed.
6. The data will appear on the display (1), percentage of fat "FAT" (pic.15), water "TBW" (pic.16), muscles mass "MUS" (pic.17) and bone mass "BONE" (pic.18).
7. The results of the calculations will be displayed (1) twice, after that the scales will switch off.

WHAT IS EXCESS WEIGHT?

It is well known that excess weight may negatively influence the health.

Excess weight may cause fatness, heart

problems, progress of atherosclerosis, osteoporosis, although the lack of weight may lead to osteoporosis too (osteoporosis is the systemic skeleton disease. Its features are: bone mass per volume decreases, bone tissue micro architectonics impairment. It issues in bones fragility, the risk of fractures because of minimal trauma or even without one.)

The physical exercises and moderate dietary regime help to get rid of excess weight. Excess weight analyzer is the ideal tool to control excess weight.

It is well known, that in case of small weight muscles and tissues contain more water than fat. And that water conducts electric current better than fat does.

The appliance can quickly measure your body fat percentage in muscles by conducting weak electric current through the body. This helps you to control your forms.

Changing of fat percentage during the day

Such factors as way of life, physical activity, health condition as well as quantity of liquid in the body influence percentage of fat.

You should take into account the following factors:

- fat percentage may be lower after taking a bath or a shower.
- fat percentage may be higher after eating or drinking.
- percentage of liquid may decrease resulting in changes of fat percentage during or after illness, as well as after physical exercises.

Water percentage

Water percentage is the amount of water in the body of a human being in percentage to total weight.

Water is of primary importance for tissues and organs of a human being, about 50-65 % of our body consist of water. You can figure out percentage of water using a formula, used by dietarians, the acquired re-

sult will be indicated as percentage of the total weight.

Thanks to the "BIA" technique used in the scales, you can control percentage of water in your body and adjust this important indicator.

Water percentage standard in female body is 45-60%, in male body - 50-65%

What is muscle mass?

Muscle mass includes weight of muscles of supporting-motor system, smooth muscles (heart and digestive system) and water, contained in these muscles.

Muscle mass should be especially controlled by people, suffering from excess weight and obesity. Weight reduction may lead to reduction of fat tissues, but only if muscle mass will be accumulated.

What is bone mass?

Human skeleton is based on bone tissues. They are hard and environmentally stable, protecting internal organs, and together with muscle tissue provide human ability to move.

The basis of bone mass measurements constitutes definition of mineral substances in the body (calcium and others).

According to the investigations, development of muscle tissues strengthens skeleton bones. This fact should be taken into account when composing dietary regime and while physical activity for bones to remain strong.

Bone mass practically does not change; nevertheless, some decrease is possible as a result of many factors (mass, height, age, gender). That is why there are no indications of the directions of changes and interpretation of measurements. There are no regulation values and recommendations.

Note: - the data in the table are of recommended character.

ENGLISH

Male	Age	% fat			% water			% muscles mass		
		Low	Medium	High	Low	Medium	High	Low	Medium	High
	20-39	5.0-7.9	8.0-19.0	19.1-45.0	38.0-54.9	55.0-70.0	70.1-85.0	45.0-72.9		82.1-95.0
	40-59	5.0-10.9	11.0-22.0	22.1-45.0	38.0-49.9	50.0-68.0	68.1-85.0	45.0-69.9	73.0-82.0	79.1-95.0
	60-79	5.0-12.9	13.0-25.0	25.1-45.0	38.0-49.9	50.0-65.0	65.1-85.0	45.0-66.9	70.0-79.0	77.1-95.0
									67.0-77.0	
Female		% fat			% water			% muscles mass		
	20-39	5.0-20.9	21.0-33.0	33.1-45.0	38.0-47.9	48.0-65.0	65.1-85.0	45.0-67.9		77.1-95.0
	40-59	5.0-22.9	23.0-35.0	35.1-45.0	38.0-44.9	45.0-65.0	65.1-85.0	45.0-64.9	68.0-77.0	75.1-95.0
	60-79	5.0-23.9	24.0-36.0	36.1-45.0	38.0-44.9	45.0-65.0	65.1-85.0	45.0-61.9	65.0-75.0	72.1-95.0
									62.0-72.0	
Male, going in for sports up to 10 hours a week		% fat			% water			% muscles mass		
	20-39	5.0-5.9	6.0-12.0	12.1-45.0	38.0-54.9	55.0-71.0	71.1-85.0	45.0-74.9		85.1-95.0
	40-59	5.0-5.9	6.0-13.0	13.1-45.0	38.0-49.9	50.0-69.0	69.1-85.0	45.0-72.9	75.0-85.0	83.1-95.0
	60-79	5.0-7.9	8.0-13.0	13.1-45.0	38.0-49.9	50.0-67.0	67.1-85.0	45.0-69.9	73.0-83.0	80.1-95.0
									70.0-80.0	
Female, going in for sports up to 10 hours a week		% fat			% water			% muscles mass		
	20-39	5.0-13.9	14.0-18.0	18.1-45.0	38.0-47.9	48.0-66.0	66.1-85.0	45.0-67.9		80.1-95.0
	40-59	5.0-15.9	16.0-20.0	20.1-45.0	38.0-44.9	45.0-65.0	65.1-85.0	45.0-65.9	68.0-80.0	78.1-95.0
	60-79	5.0-16.9	17.0-20.0	20.1-45.0	38.0-44.9	45.0-65.0	65.1-85.0	45.0-64.9	66.0-78.0	76.1-95.0
									65.0-76.0	

Time factor related to the results

Take into account that all the results of the calculations and common measuring have value only at long-term comparison of the results, short-term change of body mass within the limits of several days, as a rule, is caused by trivial loss of liquid in the tissues of your body.

Estimation of the results is made based on the changes of total weight of the body, percentage of fat, liquid in the tissues and muscles mass of the body depends on the time interval, for which these changes happened. Short-term changes of the indications during the day may differ significantly from changes of indications during a week or a month.

As a rule, short-term changes in the indications are related to the content of liquid in the tissues of the body, and medium and long-term changes are caused by changing of fat and muscles mass.

- For example, short-term decrease of weight and simultaneous increase or unchanged quantity of fat in the tissues of the body may be caused by the decrease of the liquid amount as a result of physical exercises, visiting sauna or quick loss in weight.
- Medium-term increase of weight, decrease or unchanged quantity of fat in the tissues of the body may be caused by the increase of the muscles mass.
- If you notice simultaneous decrease of total weight and fat, your dietary regime works – you lose fat in the tissues of your body.
- Ideal condition is combination of a dietary regime with different physical exercises, you can observe medium-term increase of muscles mass in the tissues of the body.

CLEAN AND CARE

- Take care using the scales. This is the precision tool, avoid temperature differences, high humidity, strikes, dust and direct sun rays.
- Keep away from heating appliances.
- Wipe the platform with a damp cloth.
- Wipe the base of the scales with a soft dry cloth.
- Do not use chemicals and abrasive agents to clean the scales.
- Do not use or keep the scales in rooms with high humidity (more than 80%), do not let it contact water or other liquids, as such a contact could affect the measurement accuracy or lead to scale malfunction.
- Keep the appliance in a dry place out of reach of children.

TECHNICAL SPECIFICATIONS

Maximal weight: 150 kg/330 lb/23 St

Scale division: 0,1 kg/ 0,2 lb/1/4lb

Batteries: 4x1,5V AAA batteries

Specifications subject to change without prior notice.

Service life of appliance - 5 years

GUARANTEE

Details regarding guarantee conditions can be obtained from the dealer from whom the appliance was purchased. The bill of sale or receipt must be produced when making any claim under the terms of this guarantee.

This product conforms to the EMC-Requirements as laid down by the Council Directive 89/336/EEC and to the Low Voltage Regulation (73/23 EEC)

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© ООО ГОЛДЕР-ЭЛЕКТРОНИКС, 2011
© GOLDER-ELECTRONICS LLC, 2011